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P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany

Instructions for Use Dentamed P&B

Light cure single component adhesive

Dentamed P&B is a simple to use all purpose light cure single component priming and bonding system which requires no mixing ("one step single bonding"). It is designed for strong bonding of composite and amalgam to enamel and dentine and of composite to nonprecious and precious metals.

Its strong adhesion to enamel or dentine proceeds on principles similar to that occurring with glassionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by carboxylic acid grafted methacrylate esters. It is therefore crucially different from other bonding systems, which achieve adhesion by means of phosphate groups, which are sensitive to hydrolysis and thus will be weakened by time. Dentamed P&B is based on acetone.

Dentamed P&B is compatible with all current brands of visible light cure methacrylate based composite restorative materials.

1. Isolation

Rubber dam is the recommended method of isolation.

2. Cavity preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of calciumhydroxide material.

4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply an etching gel (e.g. Dentamed Etching Gel) onto the enamel and dentine surfaces beginning with the enamel bevels. Leave the etching gel in place for **20 seconds**. On primary teeth and teeth high in fluoride, a **60 second** etch is recommended. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but **do not** desiccate. A slightly wet dentine surface is important for the function of Dentamed P&B. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

5. Application of Dentamed P&B

Apply Dentamed P&B generously with a brush onto the dentine and enamel surfaces for **30 seconds** with agitation. Re-apply the Dentamed P&B as needed to insure that all the surfaces are kept wet with the Dentamed P&B for the specified time. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Do not desiccate the dentine.

Note:

Do not rinse off the Dentamed P&B! If not used immediately, place dispensed Dentamed P&B in subdued light to prevent premature polymerization by incident light.

Cure the Dentamed P&B coating by exposing its entire area to a dental halogen light unit for **20 seconds** before application of a second layer of Dentamed P&B.

Note: The Dentamed P&B will **not** selfcure.

It is essential that the primed dentine and the etched enamel surfaces are dry and contaminant free for the next application.

6. Application of a second layer of Dentamed P&B

Apply again the Dentamed P&B generously with a brush onto the adhesive surfaces as described under **Point 5.** before placement of a light cure composite.

7. Restorative placement, cure and finishing

Refer to manufacturer instructions for placement, curing and finishing of light cure restorative composite.

8. Storage

Do not store above 25 °C. Do not use after expiry date.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!
For dental use only!

Vertrieb:

Dr. Ihde Dental GmbH
Erfurter Strasse 19
D-85386 Eching
Tel. +49 89 319 761-0
info@ihde-dental.de

Manufacturer:
P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany