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English:  
Instructions for use – Please read carefully

# IHDENT® Experdent® APH

Experdent® APH is a light-curing, radiopaque filling material for restorations of all cavity classes.

## Indications

Filling material for restorations of all classes of cavities.

## Contraindications

- Do not use the material if a relatively dry working area or the recommended application technology is not possible
- Do not use the material where there are known allergies to any of the components or contact allergies.
- Do not use directly on the exposed pulp.

## Basic safety information

- For dental use only.
- Keep out of the reach of children.
- Avoid contact with the skin. In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid contact with the eyes. In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

## Interactions

- Water and air containing oil can prevent the polymerization of the materials at the contact surfaces.
- The use of cationic mouthwashes, plaque disclosing agents, or chlorhexidine can lead to discoloration.
- Not be used in conjunction with materials that contain eugenol / clove oil as these inhibit polymerization and can result in staining.

## Side effects

There are no known side effects to date. Hypersensitivity to components of the material cannot be ruled out in individual cases.

## Notes for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.
- Do not remove the oxygen-inhibited layer if a further layer of filling material is to be applied, because this is required for the bond with the next layer.
- Follow the manufacturer's instructions for any other materials that are used with Experdent® APH.

## Recommended use

- Clean tooth. Remove cleaning residue with a water spray. Determine the VITA shade.
- Note:** The use of a rubber dam is recommended.
- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. Protect areas close to the pulp.
  - if necessary, etch using your preferred etching technique (e.g., Dentamed® Etching Gel, [REF](#) 5138). Observe the manufacturer's instructions for the material used.
  - Apply bonding agent in accordance with the manufacturer's instruction.
  - If necessary, apply an underfilling of low-viscosity composite or compomer. We recommend using a matrix.
  - Apply the desired shade of the Experdent® APH with a maximum layer thickness of 2 mm into the cavity. Light-cure each layer for 20 seconds.
  - Removal of excess material, contouring, and fine processing are done using diamond finishing tools and flexible polishing discs of decreasing grain size, for example. Experdent® APH can be given a high shine using polishing paste together with polishing cups or discs.

## Composition

Barium glass in a BIS-GMA matrix made up of dental resins. Pigments, additives, and catalyst. Filler content: 85 weight % = 68 vol.%. The range of variation of the inorganic filler particles is between 0.02 and 3.6 µm.

## Classification

DIN EN ISO 4049:2009, Type 1, Class 2, Group 1  
The indicated value for the depth of polymerization according to ISO 4049:2009 of Experdent® APH is 2mm. The layer thickness of 2 mm specified in the guideline is recommended.

## Storage

- Store at room temperature (15-25 °C / 59-77 °F).
- Close syringe immediately after use so that the material does not cure in the ambient light.
- Do not use after the expiration date.

## Packaging

1 syringes 3 g paste in the shades

Universal A2 [REF](#)4013

Yellow B3 [REF](#)4015

Brown A3,5 [REF](#)4017

Universal A3 [REF](#)4018

## Legends

 Read instructions

 Expiration date

 LOT Charge number

 Keep in a dry place

 Store tightly closed keep

 Temperature range from 15 °C to 25 °C

 Manufacturer

CE 0481

9999-4013-GA-01

Revision issued 2018-05-24

 **Dr. Ihde Dental GmbH**  
Erfurter Straße 19  
D-85386 Eching  
Tel. +49 (0)89 319 761 0  
Fax +49 (0)89 319 761 33  
info@ihde-dental.de  
www.ihde-dental.de

IHDENT®, Experdent® and P&B® are registered trademarks.

Vita® is a registered trademark of Vita Zahnfabrik H. Rauter GmbH.