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Manufacturer:

P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany

Instructions for Use**Ceramolar****Low shrink low abrasion light cure Nano Composite with improved handling features**

Ceramolar is a visible light cure Nano composite for anterior and posterior restorations. Beside its low shrink and low abrasion Ceramolar shows improved handling features. Ceramolar is based on urethanedimethacrylate resin and inorganic filler particles < 1.0 µm. The total filler content is 82 % by weight and 74 % by volume.

The composite meets the requirements of DIN EN ISO 4049.

Indications

Suitable for all cavity classes.

Contraindications

The placement of Ceramolar is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Ceramolar if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Ceramolar may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparation of Tooth and Colour Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and colour matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Ceramolar shade.

Ceramolar is shaded according to VITA®-Shades. It is also available in special shades and bleach shades.

Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

Best pulp protection in case of deep excavations is the photo-activated disinfection (PAD). If this is not possible, deep excavations should be covered with a thin layer of calcium hydroxide material.

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Dentamed Etching Gel and subsequent application of a light cure bonding agent (e.g. Dentamed P&B II) is recommended. Otherwise use a self etching bonding agent (e.g. Dentamed E&B) before application of the Ceramolar.

Application and Curing of Ceramolar

Apply Ceramolar in the selected shade. Ceramolar is easy to model. Place it by using a suitable instrument. Transparent matrix strips may be used.

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite (e.g. Enamelite), followed by the application of Ceramolar.

In case of small cavities Ceramolar can be applied directly.

For an optimum result apply Ceramolar in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity of at least 1000mW/cm²**, cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**

Darker shades (e.g. A3.5, B3) **30 seconds**

Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

Storage

Do not store above 25°C (78°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

APPLIC.-Tips are for single use only.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!**For dental use only!**

Vertrieb:

Dr. Ihde Dental GmbH
Erfurter Strasse 19
D-85386 Eching
Tel. +49 89 319 761-0
info@ihde-dental.de

Manufacturer:

P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany