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# Ceramolar Flow

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## Instructions for Use

# Ceramolar Flow

**Low shrink low abrasion flowable light cure nano-composite with improved mechanical properties**

Ceramolar Flow is a visible light cure high-viscous flowable nano-composite for anterior and limited posterior restorations. Beside its low shrink and low abrasion Ceramolar Flow shows improved mechanical properties.

Ceramolar Flow is based on BIS-GMA-resin and inorganic filler particles <1,0 µm.

The total filler content is 69 % by weight and 56 % by volume.

The composite meets the requirements of DIN EN ISO 4049.

### Indications

- Anterior restorations (class III, IV, V)
- Small class I and II restorations
- Lining of cavities

### Contraindications

The placement of Ceramolar Flow is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Ceramolar Flow if the patient is known to be allergic to any of the ingredients.

### Side effects

In singular cases, Ceramolar Flow may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

### Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

### Preparation of Tooth and Colour Matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and colour matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Ceramolar Flow shade.

Ceramolar Flow is shaded according to VITA®-Shades. It is also available in special shades and bleach shades.

### Cavity Preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Calciumhydroxid-Liner).

Condition and apply a bonding agent according to the manufacturer instructions. The total etch technique with Dentamed Etching Gel and subsequent application of a light cure bonding agent (e.g. Dentamed P&B II) is recommended. Otherwise use a self etching bonding agent (e.g. Dentamed E&B) before application of the Ceramolar Flow.

### Application and Curing of Ceramolar Flow

Apply Ceramolar Flow in the selected shade. Place it by using a suitable instrument. Matrix strips may be used.

For an optimum result apply **Ceramolar Flow** in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity of at least 1000mW/cm<sup>2</sup>**, cure each increment as follows:

Lighter shades (e.g. A2, A3) **20 seconds**

Darker shades (e.g. A3.5) **30 seconds**

### Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

### Storage

Do not store above 25°C (78°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

### Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

### Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

### Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

### Keep away from children!

**For dental use only!**

Vertrieb:

**Dr. Ihde Dental GmbH**  
Erfurter Strasse 19  
D-85386 Eching  
Tel. +49 89 319 761-0  
info@ihde-dental.de

Manufacturer:

P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany