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Chargennummer

**Type of preparation**

PLIWA® Roto-Pren is an autonomously acting disinfecting and cleaning agent for precision rotating instruments; it is based on bases, alcohol and amine, special solvents and material-protecting corrosion inhibitors.

**Areas of use**

PLIWA® Roto-Pren is used to clean and disinfect:

- Dental drills
- Precision rotating instruments

**Application**

Place rotating small instruments into a bur steriliser or ultrasound bath containing undiluted solution immediately after use, and without prior cleaning, in order to prevent dirt from drying onto the instruments. After removing instruments from the PLIWA® Roto-Pren rinse them carefully with water and dry them. If dirt is persistent, then brush the instruments down after removal and repeat disinfection using PLIWA® Roto-Pren. If instruments are to be sterilised following disinfection using PLIWA® Roto-Pren, then they should be rinsed and dried first. Replace used solution daily in order to achieve optimum effectiveness.

**Important:** Do not place aluminium parts in PLIWA® Roto-Pren. Do not mix PLIWA® Roto-Pren with other cleaning agents.

**Microbiology**

PLIWA® Roto-Pren is effective against:

- Bacteria
- Levurocide (Candida albicans)
- TB, mycobactericide
- BVDV / vaccinia virus Elstree strain

**Application concentrations / application time**

Bactericide high load	1 minute
Levurocide high load	1 minute
Vaccinia virus Elstree strain / BVDV	1 minute
TB, Mycobacterium terrae high load	5 minutes

**Product composition**

100 g of ready-to-use solution contains:

15 g Isopropanol / 2,97 g potassium hydroxide / 1,25 g sodium nitrate / 2,10 g N,N-bis(3-aminopropyl)dodecylamine, corrosion inhibitors

**Chemical and physical details**

Appearance of concentrate: clear fluid

pH: > 13

Density (20° C): 1,010 - 1,030 g/cm<sup>3</sup>

**Ecological information**

Biodegradability: easily biodegradable

Method: OECD 301D/EEC 84/449 C6

**Precautionary statements**

<b>H226</b>	Flammable liquid and vapour
<b>H314</b>	Causes severe skin burns and eye damage
<b>P101</b>	If medical advice is needed, have product container or label at hand.
<b>P102</b>	Keep out of reach of children.
<b>P103</b>	Read label before use.
<b>P210</b>	Keep away from heat / sparks / open flames / hot surfaces – no smoking.
<b>P243</b>	Take precautionary measures against static discharge.
<b>P260</b>	Avoid breathing dust / fume / gas / mist / vapours / spray.
<b>P280</b>	Wear protective gloves / protective clothing / eye protection / face protection.
<b>P301 / P330 / P331</b>	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
<b>P303 / P361 / P353</b>	IF ON SKIN (or hair): remove / take off immediately all contaminated clothing. Rinse skin with water / shower.
<b>P305 / P351 / P338</b>	IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.
<b>P304 / P340</b>	IF INHALED: remove victim to fresh air and keep at rest in a position comfortable for breathing.
<b>P309 / P310</b>	IF EXPOSED OR YOU FEEL UNWELL: immediately call a POISON CENTRE or doctor / physician.
<b>P403 / P235</b>	Store in a well-ventilated place. Keep cool.
<b>P501</b>	Dispose of contents / container in compliance with local / regional / national / international regulations.

**Labelling**

CE 0124 medical device Class II b / UMDNS 17-920

**Storage instructions**

Store the products at room temperature. Do not use after expiry date.

**Special note**

Always read the Instructions for Use and the Material Safety Data Sheet before use.

**Keep out of reach of children.****Delivery units**

2 L REF 10-0052

**Expert opinions:** Available on request.

Revision date 2016-05

**L.C.Pliwa®**  
 Pliwastr. 2  
 DE-34323 Malsfeld-Ostheim  
 Tel. +49 (0)5661 7317 0  
 Fax +49 (0)5661 7371 10  
 info@pliwa.de  
 www.pliwa.de

**Distribution by  
 Dr. Ihde Dental AG**  
 Dorfplatz 11  
 8737 Gommiswald  
 Tel. +41 (0)55 293 23 23  
 Fax +41 (0)55 293 23 00  
 contact@implant.com  
 www.implant.com

**Danger**

CE 0124

**Legend**

	Read instructions
	Expiration date
	LOT Charge number

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