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Chargennummer

Type of Preparation

PLIWA® Biguaton forte NEU is an instrument disinfectant based on quaternary ammonium salts, amines, biguanides and alcohol, and offering short application times and low application concentrations. This gentle disinfection and cleaning agent cleans over a wide degree of efficiency. PLIWA® Biguaton forte NEU contains corrosion inhibitors that protect instruments. Because it is simple to use, it enables cleaning and disinfection in one single process.

Areas of use

PLIWA® Biguaton forte NEU is a highly effective instrument disinfectant for dental and medical use. Suitable for medical instruments, laboratory equipment made from glass, porcelain, metal, rubber, plastic and so on.

Application

PLIWA® Biguaton forte NEU is supplied as a concentrate and should be diluted using water at 30–40° C, then kept at the required concentration. Pre-clean instruments immediately after use, and place them into the stock solution. Any surfaces and hollow spaces which are to be disinfected must be wetted completely by the stock solution. After the disinfection period, rinse the instruments thoroughly and dry them. In order to retain the disinfectant's microbiological efficacy, please do not exceed the application times or concentrations. Suitable for all circulation processes and ultrasound devices. We recommend that you replace the solution daily in order to achieve optimum effectiveness.

Application concentrations / application time**Application time as defined by Association of****Applied Hygiene:**

	1%	5 minutes
Bactericide, levurocide (C. albicans) low load	1% 0,5% 0,25%	5 minutes 15 minutes 30 + 60 minutes

tuberculocide (M. terrae) low load	3% 2% 1%	15 minutes 30 minutes 60 minutes
---------------------------------------	----------------	----------------------------------------

Virucide

limited viricide HCV (testvirus BVDV) 0,5% 5 minutes

HBV/HIV/H1N1/H5N1
HSV incl. SARS-CoV
Vaccinia elstree strain 1% 15 minutes

Viricide
MVA modified Vacciniacirus Ankara 0,25% 15 min clean
1% 15 min dirty
conditions

Product composition

100 g PLIWA® Biguaton forte NEU contains:

N,N-Didecyl-N-methyl-poly(oxyethyl)-ammonium propionate, Propandiamine, N-(3-aminopropyl)-N-dodecyl, Coco propylenediamine-1,5-bis-guanidinium acetate, Fatty alkylamine

Chemical-physical details

Appearance of concentrate: dark blue fluid
pH (concentrate): approx. 10,5 (± 0.5)
pH (1% solution): 9,4
Density (20° C): 0,940–0,999

Ecological information

Biodegradability: easily biodegradable
Method: OECD 301D / EEC 84/449 C6

Interactions

If PLIWA® Biguaton forte NEU comes into contact with other stock solutions or preparations, it may become less effective and may discolour. Important: Instrument baths must be cleaned thoroughly prior to first-time use of PLIWA® Biguaton forte NEU. They must not contain any remnants of other instrument solutions.

Precautionary statements

H302	Harmful if swallowed
H314	Causes severe skin burns and eye damage
H400	Very toxic to aquatic life
P260	Do not inhale dust / fume / gas / mist / vapours / spray.
P280	Wear protective gloves / protective clothing / eye protection / face protection.
P273	Avoid release into the environment.

P301 / P330 / P331

IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303 / P361 / P353

IF ON SKIN (or hair): remove / take off immediately all contaminated clothing. Rinse skin with water / shower.

P305 / P351 / P338

IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.

P301 / P312

IF SWALLOWED: call a POISON CENTRE or doctor / physician if you feel unwell.

P304 / P340

IF INHALED: remove victim to fresh air and keep at rest in a position comfortable for breathing.

P309 / P310

IF EXPOSED OR YOU FEEL UNWELL: immediately call a POISON CENTRE or doctor / physician.

P501

Dispose of contents / container in compliance with local / regional / national / international regulations.

Labelling

CE 0124 Class II b / UMDNS 17-920

UN 1903

Storage instructions

Do not use after expiry date. In extremely cold temperatures, cloudiness and precipitation can appear in instrument disinfection solutions that have been mixed. This can be dispersed easily by heating to above 30°C. Cloudiness does not mean a loss of effectiveness. Keep container closed tightly. Protect against heat, direct sunlight and frost.

Special note

Always read the Instructions for Use and the Material Safety Data Sheet before use.

Keep out of reach of children.**Delivery units**

2 L REF 10-0002

Expert opinions: Available on request.

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Dr. Ihde Dental AG