Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Enamelite®

Enamelite® is a low-viscosity light-curing microhybrid composite for anterior and small posterior restorations and for pit and fissure sealing.

Essential components

Enamelite® is based on methacrylates and inorganic fillers with particle sizes between 0.05 and 1 µm. The total filler content of the material is 61 % by weight or 41 % by volume.

Application areas

Black class III, IV, V cavities, Sealing

Contraindications

None reported. Susceptible persons may exhibit hypersensitivity reactions to Enamelite®. If this occurs, discontinue the use of the product immediately. For affected patients, the Enamelite® filling should be removed using rubber dam protection in order to prevent the inoculation of the mucosa with filling residue.

Adverse reactions

Direct contact with the dental pulp may cause irritation. In deep cavities, the dental pulp must therefore be protected by lining the cavity floor with a thin layer of calcium hydroxide liner (REF 8768).

Side effects

Do not use in connection with preparations containing eugenol. Eugenol may inhibit the polymerization process and can result in discolouration of the filling. Storing the product near products containing eugenol may also be detrimental unless the latter are tightly sealed. For this reason, never expose unpolymerized composite material to products containing eugenol.

Application

Preparing the tooth

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Enamelite® shade.

Polymerization

20 seconds of light polymerization are recommended for the A2 shade and 30 seconds for the A3 and A3,5 shades.

Note

To minimize polymerization shrinkage, polymerize the material in 1 mm increments for all restorations whose depth exceeds 1 mm.

Warning

The consistency of the composite material may not be modified by adding bonding. The composite material must not be diluted with resins. Avoid contact of the unpolymerized Enamelite® with the unprotected skin and eyes.

Keep out of the reach of children.

Intended for use by dentists only.

Storage / Shelf life

Do not store above 25 °C. Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Package sizes

Syringe 5 g	A2	REF 4050
Shade: universal	A3	REF 4051
Shade: brown	A3,5	REF 4054
3 application tips		

Legends

	Read instructions
	Expiration date
	LOT Charge number
	Store protected against the sun
	Keep in a dry place
	Store tightly closed keep
	Temperatures up to 25 °C
	Manufacturer

CE 0029

DIN EN ISO 4049

9999-4050-03

Revision issued 2015-07

Dr. Ihde Dental AG

Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH

Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT® and Enamelite® are registered trademarks.

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