Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Priginal® H

Priginal® H is a modern, self-curing two component hybrid composite for dental fillings based on chemically bound filling particles.

Constituents

Priginal® H is based on Bis-GMA-resin and anorganic filling particles in the range of 0.05 - 5 micron. Portion of the total filling material: 61 % volume of the anorganic filling material: 44 %.

Indications

- Fillings of the cavity classes I and II
- Aesthetic correction in visible area
- Correction of morphological anomalies in enamel
- Repair of plastic and ceramic veneers
- Substructures under crowns and bridges

Use

The use of the acid etching technique allows a substance saving preparation of the tooth. Prior to the preparation, cleanse the tooth with fat-free pumice paste and water. Using a cofferdam eases the drying. Prepare the cavity with a possibly low amount of tooth substance loss. The enamel surface surrounding the cavity is chamfered. These enamel areals are etched with IHDENT® Etching Gel (REF 5138) for 20 - 30 seconds. Please note the instructions of the etching substance! Afterwards the IHDENT® Etching Gel is washed off with a pure water jet. The tooth surface is dried and should have a dull shine. Contamination with saliva should now be avoided. IHDENT® P&B® is mixed in a 1:1 ratio and applied to the etched cavity area while avoiding contact with the tooth, and afterwards blown out to a thin layer with the air blower. The base and active paste are mixed in a 1:1 ratio. Plastic spatulas are preferably used for this since otherwise debris from the metal spatula can enter the composite due to the larger filling particles. Mix the paste thoroughly for 40 seconds. Correctly mixed material should be free of reams and appear homogenous. The application of the composite can be performed with metal instruments. The mixture is processable approximately 1:30 minutes. The forming occurs with plastic matrices or crowns. If no forming aids are used, the top layer of the composite cannot cure. The unpolymerized layer must in any case be removed by polishing. The curing occurs within 7 minutes after beginning mixing. It may not be disturbed.

Polishing / Curing

Priginal® H fillings can be easily polished. After the rough contouring with finishing diamonds or stones, the polishing preferably occurs with polishing discs of interchangeable roughness. Sufficient cooling is to be ensured. Priginal® H can be applied in several layers. An intermediate polishing is thereby not necessary.

The following factors can influence the stability of Priginal® H:

- Contamination of the pastes among themselves as well as with fat, dust etc.
- Heat impact, also brief.

Side effects

Especially if even a thin layer of dentin remains to the vital pulpa, a toxic effect of the composite can not be excluded. Therefore an underfilling is always indicated in order to protect the vital pulpa. Materials containing eugenol may not be used for this purpose.

Plastic monomeres can cause allergies in sensitive persons. Therefore, sufficient curing is to be ensured. In case an allergy to a substance in the composition is known, the use is contraindicated. Should the composite fillings be removed upon occurrence of an allergy, the removal may only be performed with the use of a cofferdam in order to avoid inoculation of the plastic particles.

A prophylactic removal without cofferdam is strictly contraindicated. Contact of the unprotected skin surface with unpolymerized composite is contraindicated in order to avoid allergy formation.

Contraindications

None known.

Interactions

Phenolic substances impair the polymerization. Therefore, underfilling materials containing eugenol may not be used.

Shelf life / Storage

Store at room temperature (15 - 25 °C / 59 - 77 °F). Do not use after expiration date. Expiration date and production no.: see barcode.

Precautionary statements

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

Keep out of reach of children. Only to be used by the dentist.

Package size

14g Base paste
14g Active paste (color A2)
Mixing spatula **REF 4021**

Recommended bonding:

IHDENT® P&B® **REF 5120**

Legends

	Read instructions
	Expiration date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperature range from 15 °C to 25 °C
	Manufacturer

CE 0029

ISO 4049

Warning

- H315**
Causes skin irritation.
- H317**
May cause an allergic skin reaction.
- H319**
Causes serious eye irritation.
- H335**
May cause respiratory irritation.

9999-4021-03

Revision issued 2015-08

Dr. Ihde Dental AG
Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH
Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT®, Priginal® and P&B® are registered trademarks.

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