Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Triomolar®

Light-curing micro-hybrid composite for anterior and posterior restorations.

Triomolar® is classified as a light-curing microhybrid composite. It is characterized by the following outstanding properties:

- Triomolar® is extremely abrasion-resistant and offers excellent edge stability and strength.
- Triomolar® can be polished to permanent high gloss.
- Thanks to its optical properties, properly processed Triomolar® is nearly indistinguishable from natural tooth substance.
- Triomolar® allows for excellent long-term restorative results even in areas with high loads.
- Triomolar® is strongly radiopaque.

Essential components

BIS-GMA resins and inorganic fillers having a particle size of 0,05 to 0,9 µm. The total filler content is 81% (weight percent) and 63% (volume percent).

Indication

- Class I–V fillings where a proper indication exists.
- Correction of morphological anomalies in enamel and dentine regions.
- Repair of resin or ceramic veneers.
- Crown and bridge substructures.

Adverse reactions

Toxic effects of composites on the dental pulp cannot be ruled out, particularly if only a thin layer of dentine remains between them and the vital pulp. To protect the vital dental pulp, a cavity liner containing e.g. calcium hydroxide, REF 8768 is indicated – particularly if the residual dentine is less than 1 mm thick. In these cases, cavity liners containing eugenol are contraindicated.

Contraindications

In sensitive persons, monomers of composites can induce hypersensitivity (allergy). As soon as the first symptoms appear, the use of the substance must be stopped. Triomolar® must not be used in patients with a known history of allergic reactions to one or more of the substances contained in Triomolar®. If composite fillings in hypersensitive persons must be removed, this should be done using rubber dam isolation to prevent the inoculation of particles into the mucosa. Any removal without rubber dam protection is strictly contraindicated. To avoid development of an allergy, contact of the substance with the unprotected skin is contraindicated.

Interactions

Phenolic substances prevent polymerization. Cavity liners containing eugenol should therefore be avoided. Eugenol can inhibit the composite curing process and may result in colour changes. Even storing open containers with products containing eugenol near composite materials may result in damage. Do not permit non-polymerized composite materials to be exposed to products containing eugenol. Do not use any resin to adjust viscosity of composite restorative material.

Application

Preparation of tooth

Clean the tooth with a fluoride-free polishing paste (e.g. fluor of pumice with water) prior to preparation and colour matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Triomolar shade.

Cavity preparation

After isolation (best with the rubber dam) prepare the cavity with minimal tooth reduction. Best pulp protection in case of deep excavations is the photoactivated disinfection (PAD). If this is not possible, deep excavations should be covered with a thin layer of hard setting calcium hydroxide material.

Condition and Application

The acid-etching technique and the use of bonding are recommended because this allows conservative preparation, which preserves tooth substance.

The following bonding systems are recommended:

Dentine and enamel: IHDENT® P&B® (REF 5120); single Component primer and bond

In principle, Triomolar® can be used with all modern bonding systems.

Preparation of dentine and enamel

Pre-treat the dentine and enamel with 37% phosphoric acid for 15 seconds. Rinse the tooth with water and dry the tooth surface. Highly fluoridated enamel areas can be etched with phosphoric acid for up to 70 seconds to improve the retentive etching pattern. If the residual thickness of the dentine is low (particularly if under 1 mm), place a cavity liner. Etching for longer than 15 seconds in the dentine region is not recommended under any circumstances.

Apply a suitable bonding agent to the rinsed and dried tooth surface. Light-polymerization of the bond layer is always recommended for the resin tags penetrating the dentinal canals to harden completely.

Application and Curing of Triomolar®

Apply Triomolar® in the selected shade by using a non-metallic instrument. Contour and polymerize using a VLC lamp as usually found in a dental office. Triomolar® is easy to model. Transparent matrix strips may be used. It is recommended to build up the restoration in layers in the molar region. In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement or a light cure flowable composite, followed by the application of Triomolar®. In case of small cavities Triomolar can be applied directly. For an optimum result apply Triomolar® in layers of max. 2 mm. Light cure each increment separately.

By using a polymerization unit with a light intensity of at least 1000 W/cm², cure each increment as followed:

Lighter shades (such as A1, A2, B1) **20 seconds**

Darker shades (such as A3.5 or B3) **30 seconds**

To minimize the effect of the polymerization shrinkage, it is recommended to polymerize in layers if the restoration is deeper (higher) than 2.5 mm.

Finishing / Polishing

Triomolar® can be polished to a perfect high-gloss finish. Use either polishing discs or rubber polishers. To protect the periodontal tissues, remove excess material in the proximal area with polishing strips.

Storage / Shelf life


Stability in storage: For the expiry date and batch number, see the product label. Storage temperature < 77 °F (< 25 °C). Lower storage temperatures improve storage stability. Store Triomolar® protected from light, in a dry place.

Keep this and all medicines out of the reach of children. To be used by dentists only.

Package sizes

Colours	Syringe, 3.5 g	20 Tips, 0.3 g
	REF	REF
Universal A2	4003	4033
Light B1	4004	4034
Yellow B3	4005	4035
Grey C3	4006	4036
Brown A3.5	4007	4037
Universal A3	4008	4038

Legend

-  Read instructions
-  Expiration date
-  LOT Charge number
-  Store protected against the sun
-  Keep in a dry place
-  Store tightly closed keep
-  Temperatures up to 20 °C
-  Manufacturer

CE 0029

DIN EN ISO 4049

9999-4003-03

Revision issued 2015-08

Dr. Ihde Dental AG
Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH
Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT®, Triomolar® and P&B® are registered trademarks.

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