Deutsch:

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English:

Instructions for use – Please read carefully

## IHDENT® CemBase® ivory (Powder)

CemBase® is a water-mixable glass-ionomer cement used for cavity liners, core build-ups and cementation in prosthetic, conservative and orthodontic dentistry. It is quick setting and fluoride releasing.

### Essential components

Aluminium silicate fluoride, polyacrylic acid, iron oxide pigments (known)

### Areas of application

Cementation of prosthetic restorations and orthodontic appliances (e.g. crowns, bridges, inlays / onlays, bands) Cavity liners

### Interactions

None known

### Side-effects

Severe irritation of the pulp may occur.

**Note:** It cannot be ruled out that CemBase® may cause sensitisation in sensitive patients. CemBase® should no longer be used if these symptoms occur. CemBase® should not be used if there is a suspected allergy. If it is necessary to remove the CemBase® in such cases, a rubber dam should be placed before removal to prevent components of the material infiltrating the oral mucosa. It is essential to apply calcium hydroxide (REF 8768) with deep cavities and minimal residual dentine layer (dentine third near the pulp) to prevent possible pulp reaction.

Over dispensing the water may lead to increased solubility of the cement and consequently pulp reaction. Hydrostatic pressure in close fitting crowns can also cause pulp reaction. Excessive over filling of crowns should therefore be avoided.

### Contraindications

Do not apply directly to exposed pulp. Use is only indicated if it is possible to ensure relative isolation.

### Caution

Avoid contact with the eyes.

### Application - Preparation

Carefully clean and dry the enamel, dentine and metal surfaces to ensure an optimal bond of CemBase®. Avoid overdrying, e.g. with an air syringe. It is sufficient to dry the surface using a pledget. A moisture film, as is usual with vital teeth, does not impair setting of the cement. Condition the dentine surface for 10 seconds using Ihdent® GIZ® Conditioner (REF 2330) to improve the bond strength. Then rinse off the conditioner and dry the cavity. Avoid recontamination with water or saliva after this cleaning stage.

### Pulp protection

It is essential to apply calcium hydroxide (REF 8768) with deep cavities and minimal residual dentine layer (dentine third near the pulp). The use of a conditioner is contraindicated in this area.

### Dispensing

Minimum 2 scoops of powder and 1 drop of water. The powder and water should be mixed on a mixing pad or clean glass plate. To ensure accurate dispensing, fill the spoon liberally with powder and wipe off the excess on the neck of the bottle when removing the spoon from the bottle. Normal tap water can also be used. Hold the water bottle vertically when dispensing. Adhere strictly to the powder / water ratio for cementation. If a lower viscosity is required, e.g. for cavity liners, this can be attained by increasing the proportion of water accordingly. Mix on a cooled glass plate to compensate for the resulting shorter working time.

### Mixing

Abrasion-resistant metal cement spatulas can be used for mixing. It is not necessary to use plastic spatulas. The powder and water should be mixed on a mixing pad or glass plate. The powder should be worked into the liquid to produce a viscous consistency. (The required consistency for cementation is attained when the cement only just drops from the spatula.) We recommend mixing portions of powder into the liquid at a time (e.g. in 3 portions) to attain better wetting.

### Working time

The following times apply for the individual working stages at a room temperature of 23 °C and 50 % humidity:

Mixing time **15 seconds**  
Working time **2:25 minutes**  
Setting time **2:25 minutes**  
Total setting time **4:45 minutes**

Higher temperatures shorten the working time and lower temperatures extend the working time. Exceeding the working time will result in a loss of adhesion to the enamel and dentine.

### Application

Apply a thin layer uniformly to the prepared tooth and crown margins and fit the restoration to its final position without applying pressure. Avoid excessive overfilling. Excess cement should only be removed once the cement has fully set. Protect the filling against water penetration with IHDENT® GIZ® Varnish (REF 2440). These cement segments should be only be removed once the cement has fully set.

**Note:** CemBase® adheres to metal instruments and should therefore be washed off using cold water before setting.

### Storage

The maximum storage temperature is 25 °C. Store dry, tightly close the container after use. Refer to the barcode for batch number and expiry date.

**Keep out of reach of children. Intended only for use by dentists.**

### Pack sizes

Pack with 50 g powder, accessories **REF 8520**

### Legend

	Read instructions
	Expiry date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperatures up to 25 °C
	Manufacturer

CE 0029

ISO 9917:2007

9999-8520-03

Revision issued 2015-08

**Dr. Ihde Dental AG**  
Dorfplatz 11  
CH-8737 Gommiswald/SG  
Tel. +41 (0)55 293 23 23  
Fax +41 (0)55 293 23 00  
contact@implant.com  
www.implant.com

**Distribution by**

**Dr. Ihde Dental GmbH**  
Erfurter Str. 19  
D-85386 Eching/München  
Tel. +49 (0)89 319 761-0  
Fax +49 (0)89 319 761-33  
info@ihde-dental.de  
www.ihde-dental.de

IHDENT®, CemBase® and GIZ® are registered trademarks.

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