IHDE[®] Exakt KN 700 / IHDE[®] Exakt S 1000 / IHDE[®] Genau 2000 / IHDE[®] Universal Pastenhärter

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Elastomeric precision impression material on polysiloxane basis, condensation curing.

Product description

IHDENT® Exakt KN 700, IHDENT® Exakt S 1000, IHDENT® Genau 2000 and Universal paste hardener are precise dental impression materials based on C silicone basis with outstanding cost effectiveness to meet the highest demands. The products are produced from selected raw materials using the most modern manufacturing technology. Special fillers in the putty material make mixing easy and provide a perfectly homogeneous result. All correction materials offer consistency, hardness, thixotropism and hydrophilia to coordinate with the various applications.

Indication

Crowns and bridges, partial crowns, inlays, onlays, total and partial dentures.

Impression technique

IHDENT® Exakt KN 700:

Preliminary impression material for simultaneous putty wash technique (1-step) and putty wash technique (two-step).

IHDENT® Genau 2000:

Functional impressions, rebasing, preliminary impression material for simultaneous putty wash technique (1-step).

IHDENT® Exakt S 1000:

Wash material for putty wash technique and simultaneous putty wash technique.

Dispensing / Processing

IHDENT® Exakt KN 700

Dispensing with Universal paste hardener:

Remove the required amount of IHDENT® Exakt KN 700 with the measuring spoon, remove the excess and place on the mixing block. Dispense 2.5 cm activator paste per measuring spoon of putty (equivalent to 2.5 divisions on the scale on the mixing block, see dispensing table) and add to the putty. Pull the putty out to form a small cord, fold it together again and knead with the tips of the fingers. Repeat this procedure several times until a homogeneous shade is obtained, but for at least 30 seconds. Observe the mixing time. Close containers and tubes immediately after use.

IHDENT® Genau 2000 / IHDENT® Exakt S 1000

Dispensing with Universal paste hardener:

Express the required amount of IHDENT® Genau 2000 or IHDENT® Exakt S 1000 in the form of a cord on the mixing block. Then place half that length of activator paste beside it (see dispensing table). Pick up the activator with a spatula, spread it evenly over the other component and mix with a spatula for at least 30 seconds until a homogeneous colour mix is obtained. Smooth the mixed material with the spatula to remove any air bubbles.

Over- and underdispensing

The total curing time can be affected by adding more or less of the activator. In general, a greater amount produces faster curing and a smaller amount means slower curing. Over- or underdispensing the activator by up to 25 % is possible.

Impression procedure

When using a suitable tray adhesive, please follow the manufacturer's instructions for use. Place the tray loaded with the silicone material into the mouth, press it on briefly and then hold it in situ until it has set fully. To check the degree of hardening always check the material in the mouth. When using the putty wash technique (two-step) the first impression must always be cleaned and dried in order to ensure curing with the putty wash material.

Disinfection

The impressions can be disinfected with 2 % glutaraldehyde solution. After removing from the mouth, rinse the impression under running water for 15 seconds. Then disinfect the impression by dipping it completely into disinfectant. Always follow the disinfectant manufacturer's instructions.

Electroplating

The impressions can be electroplated with copper or silver in the usual baths.

Model fabrication

Filling the mould is possible after 30 minutes and up to 7 days. There are no other time restrictions. All class 3 and 4 gypsums that comply with ISO 6873 can be used. The impressions can be cast several times.

Working and safety instructions

The above-listed products must be used only in accordance with their instructions for use. Any other use, which is not in agreement with these instructions for use, is the sole responsibility of the user.

- Close containers and tubes immediately after use.
- Do not leave any impression material residues in the oral cavity.
- Contact with clothing should be avoided as cured silicones are chemically stable and produce marks that cannot be removed.
- To prevent incompatibilities, do not combine products with materials from other manufacturers.
- Do not use IHDENT® Exakt KN 700, IHDENT® Exakt S 1000, IHDENT® Genau 2000 and Universal paste hardener together with polyether, polysulphides or addition-curing silicones.
- Storage: Store between 5 °C and 27 °C.
- Do not use after the expiry date.
- Processing temperature: Room temperature between 18 °C and 25 °C.

Warnings for Universal paste hardener

- **Inhalation:** Bring the person into the fresh air and treat with oxygen if necessary.
- **Eye contact:** Avoid eye contact with the activator or impression material, if necessary, irrigate the eyes immediately with plenty of clear water and consult an ophthalmologist promptly.
- **Skin contact:** Avoid prolonged skin contact with the activator. After prolonged contact wash the skin with soap and water. We recommend to wear suitable protecting clothing and gloves.
- **Swallowing:** After swallowing, rinse out the mouth and drink plenty of water. Always contact a doctor.

- **Allergic reactions:** Allergic reactions can occur in sensitive person. In case of doubt, consult a dermatologist or allergologist. If acute allergic reactions should occur during the treatment, treatment with the product should be stopped immediately.

Precautionary statements for Universal paste hardener

- P101** If medical advice is needed, have product container or label at hand.
- P102** Keep out of reach of children.
- P103** Read label before use.
- P261** Avoid breathing dust / fume / gas / mist / vapours / spray.
- P280** Wear protective gloves / protective clothing / eye protection / face protection.
- P281** Use personal protective equipment as required.
- P321** Specific treatment (see on this label).
- P405** Store locked up.
- P501** Dispose of contents / container in accordance with local / regional / national / international regulations.

Hazard statements for Universal paste hardener

- H315** Causes skin irritation.
- H317** May cause an allergic skin reaction.
- H361** Suspected of damaging fertility or the unborn child.
- H412** Harmful to aquatic life with long lasting effects.

Warning

Please note the warning in the safety data sheet also.
Information at: www.implant.com
Errors and changes reserved.

Package sizes

IHDENT® Exakt KN 700

Tin 900 ml = 1.6 kg
Bucket 10 kg

REF 7100
REF 7110

IHDENT® Exakt S 1000

Tube 150 ml

REF 7200

IHDENT® Genau 2000

Tube 150 ml

REF 7300

Universal paste hardener

Paste hardener blue, tube 60 ml
Paste hardener red, tube 60 ml

REF 7400
REF 7410

Legends

Read instructions

Expiration date

LOT Charge number

Store protected against the sun

Keep in a dry place

Store tightly closed keep

Temperature range from
5 °C to 27 °C

Manufacturer

9999-7100-7200-7300-7400-03

Revision issued 2015-08

IHDENT® is a registered trademark.

Technical data and product characteristics

Table 1:

Technical data	IHDENT® Exakt KN 700	IHDENT® Genau 2000	IHDENT® Exakt S 1000
Consistency / EN ISO 4823	Kneadable / Type 0	Heavy body / Type 1	Medium body / Type 2
Colour	White	Green	Blue
Mixing time	30 sec.	30 sec.	30 sec.
Working time incl. mixing time*	1 min. 45 sec.	1 min. 45 sec.	1 min. 45 sec.
Time in mouth	4 min.	4 min.	4 min.
Total setting time	5 min. 45 sec.	5 min. 45 sec.	5 min. 45 sec.
Deformation under pressure	2 %	5 %	6 %
Recovery from deformation	98,8 %	99,5 %	99,7 %
Linear dimensional change	0,2 %	ca. 1 %	ca. 1 %
Flavour	mint	mint	mint

* At 23 °C / 73.4 °F. The time specifications in the upper table relate to a relative humidity of 50 +/- 10 %. In general, higher temperatures accelerate while lower temperatures delay the setting.

Dispensing table

Table 2:

	Universal paste hardener
IHDENT® Exakt KN 700  1 $\triangle$ 12 ml (22,4 g)	2,5*  $\triangle$ 0,53 ml (0,5 g)
IHDENT® Genau 2000  2* $\triangle$ 1,22 ml (1,7 g)	1*  $\triangle$ 0,21 ml (0,2 g)
IHDENT® Exakt S 1000  2* $\triangle$ 1,16 ml (1,5 g)	1*  $\triangle$ 0,21 ml (0,2 g)

* 1 mixing block measurement unit: 1 cm (0.39 inch)

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Dr. Ihde Dental AG