Deutsch:

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English:

Instructions for use – Please read carefully

**IHDENT<sup>®</sup> Ihdent alloy<sup>®</sup> K**

Ihdentalloy<sup>®</sup> K has passed the cytotoxicity test according to the norm ISO 10993-5. Ihdent alloy<sup>®</sup> K is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, down to two tenth of millimetre. Its molecular structure allows you to obtain smooth, compact surfaces with little oxide formation even after repeated oven treatment. Its thermal expansion coefficient is ideal for every ceramic of the latest generation.

Ihdentalloy<sup>®</sup> K is produced in conformity to standards ISO 9693:1999 and ISO 22674:2006, it is highly corrosion resistant and it is free of toxic elements like beryllium, indium and gallium.

**Essential components**

Co:	64 %
Cr:	21 %
Mo:	6 %
W:	6 %
Others	Si, Fe, Mn

**Moulding**

To avert shrinkage blisters and cracks avoid hollows, sudden elbows and impurities. Wax walls should be 0,5 mm thick, so that after manufacturing process metal wall thickness is at least 0,3 mm. Transitions from wax to metal should be realised without large differences in thickness of the modellation.

**Casting Sprues**

Bridges: length of the strut should match the bridge body Ø:

- Connections between cone and beams: Ø 3,5 to 4 mm
- Beams: Ø 3,5 to 4 mm
- Connections between beams and crown: Ø 2,5 mm; min. length 3,5 mm

Single crowns: they can be directly sprued, casting sprue depends on object dimensions; we suggest round casting sprues: Ø 3,0 mm to 4,0 mm.

**Casting**

Melt Ihdent alloy<sup>®</sup> K in pre-heated crucibles used for this alloy only. Do not glaze and do not apply investing technologies under pressure. Reach a temperature of 850-950 °C; maintenance time in pre-heated oven depends on investment characteristics. Check that objects to be casted are positioned outside the center of the muffle.

Alloy quantity depends on: wax framework weight with sprues x metal specific weight.

- High-frequency melting equipments: when last ingots have melted in the pre-heated ceramic crucible and incandescence shadow has disappeared, start melting process immediately. Do not use any flux.
- Open-flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux.

Sample values for flame regulation:

- Acetylene 0,4 bar / Oxygen 2 bar
- Propane 0,2 bar / Oxygen 2 bar
- Leading pressure methane / Oxygen 2 bar

Overheating the melted material could cause shrinking hollows, microporosity and building of a coarse-grained structure, and cause bridges to break or ceramic investments to crack.

**Manufacturing**

Let the framework cool down at room temperature, remove the investment and sandblast with aluminium oxide 100 to 250 µm at a pressure of 3 - 4 bar.

Continue manufacturing process using fraises for hard metal. The worked framework has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3 - 4 bar and then steameclean. After cleaning, the framework should not be touched anymore and should be held with haemostatic forceps only.

**Oxide firing**

Fire for 5 min. under vacuum at 950 - 980 °C (10 °C more than the ceramics base material firing temperature). After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminium equipment 110 to 150 µm at a pressure of 2,5 - 3 bar, then steam. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.

**Ceramic investment**

Do not use any Bonder.

It is possible to use all ceramics normally available on the market for CoCr alloys. Apply manufacturing instructions of the ceramic producer.

**Re-use of sprues and cones**

The best results are obtained with pure Ihdent alloy<sup>®</sup> K. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

**Brazing and welding**

For brazing and welding it is recommended to use a product adjusted to the ingredients of Ihdent alloy<sup>®</sup> K.

**Waste disposal**

The processing scrap must be disposed of as special waste in accordance with the EC directives 91/156/CEE on waste, 91/689/CEE on hazardous waste and 94/62/CE on packaging and packaging waste and in compliance with national legislation in force on the subject.

**Safety instructions**

1. Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. An adequate test is thus advisable.
2. Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
3. Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and sandblasting.
4. We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
5. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results.

**Important note**

Use before: not relevant. To determine the batch number, check the barcode on the packaging.

**Keep away from children! For dental use only!****Delivery units**

Package of 100 g

REF 7620

**Legends**

	Read instructions
	Expiration date
	LOT Charge number
	Manufacturer

CE 0029

**ISO 22674:2006**

9999-7620-03

Revision issued 2015-07

**Dr. Ihde Dental AG**  
Dorfplatz 11  
CH-8737 Gommiswald/SG  
Tel. +41 (0)55 293 23 23  
Fax +41 (0)55 293 23 00  
contact@implant.com  
www.implant.com

**Distribution by**

**Dr. Ihde Dental GmbH**  
Erfurter Str. 19  
D-85386 Eching/München  
Tel. +49 (0)89 319 761-0  
Fax +49 (0)89 319 761-33  
info@ihde-dental.de  
www.ihde-dental.de

IHDENT<sup>®</sup> and Ihdent alloy<sup>®</sup> are registered trademarks.

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