Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Ihdentalloy® MG

Ihdentalloy® MG is a top-quality cobalt-chromium-molybdenum casting alloy. It is free of nickel and beryllium and exhibits excellent casting fit and flowability and strength for thin-walled casting high objects (down to 0.2 mm). The molecular structure facilitates a smooth, firm polish with minimal oxide formation.

Essential components

Co: 62 %
Cr: 29 %
Mo: 8,5 %
Others: C, Si, Fe, Mn

Application areas

Cast frameworks, combination restorations, crowns and bridges

Application

Pre-heating

Preheat the muffle as instructed by the manufacturer of the investing compound. Before casting, set the furnace temperature for between 950 °C and 1050 °C.

Melting and casting using the induction centrifuge

Place a few metal cubes (corresponding to the size of the casting object) in the crucible and set the unit for maximum heating power to accelerate the melting process. Start the centrifuge as soon as the metal has melted and the outer oxide layer shows signs of tearing.

Melting and casting using an open flame and centrifuge

Fire up the burner and adjust the flame so that its core is blue and about 4 mm in length. Place the metal cubes in the pre-heated crucible and evenly heat the metal, applying the flame with circular movements. Start the centrifuging process as soon as the metal has melted and the outer oxide layer shows signs of tearing. Sandblast the framework with 120-150 µm aluminium oxide. Use the standard technique for subsequent processing.

Physical and mechanical properties

Melting range: 1177 °C - 1391 °C
Density: 8.3 g/cm³
Vickers hardness: 396 HV10
Yield strength: 645 MPa
Modulus of elasticity: 210 GPa
Expansion: 4.7 %

Re-casting sprues

For best results, use only pure Ihdentalloy® MG. The sprues can be re-used together with an equal amount of virgin alloy, provided these are both from the same batch.

Brazing and welding

For brazing and welding it is recommended to use a product adjusted to the ingredients of Ihdentalloy® MG.

Safety rules

The processing scrap must be disposed of as special waste in accordance with the EC directives 91/156/EEC on waste, 91/689/EEC on hazardous waste and 94/62/EC on packaging and packaging waste and in compliance with national legislation in force on the subject.

Warning notes

- Cobalt-chrome based alloys may seldom cause dermatitis on sensitive subjects. We advise doing a contact test (i.e. Patch test).

- Verify before the application of the prosthesis if there are other implants in the oral cavity of the patient. If different metals are present a "pile effect" may occur.
- Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.
- Metal dusts and smoke are dangerous for the health. During casting and polishing use an adequate suction system.
- It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results.

Important note

Use before: not relevant. To determine the batch number, check the barcode on the packaging.

Keep away from children! For dental use only!

Delivery units

Package of 1'000 g REF 7625

Legends

-  Read instructions
-  Expiration date
-  LOT Charge number
-  Manufacturer

CE 0029

ISO 22674:2006

9999-7625-03

Revision issued 2015-07

Dr. Ihde Dental AG

Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH

Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT® and Ihdentalloy® are registered trademarks.

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